

# ECO STAND UP SLING



- Only available in jersey
- One size fits all



REF:

SXPECO2-PF-TU



The sling is positioned at waist level on the WAYUP and the adjustment is made with the sling with the patient's arms positioned at 120°.

The sling is positioned at waist level and/or sub-axillary on the EASYLEV according to the patient's tone.

## HOOK ADJUSTMENT



Two levels of hooks

## SLINGS POSITIONNING



UNCROSSED

Decent bust tonicity



CROSSED

Lack of bust tonicity