




Wall Mounts

Independent communication for each stage of the day – this is our mission. We strive for the most versatile mounting in any situation. Often a wall mount is the perfect amendment to wheelchair, table and floor mounts. It can also be simply used as a parking for your mounts.



M3D

Rehadapt wall mounts with base plates for mounting Monty 3D tubes. Perpendicular or parallel to the wall - with 3D arm as a complete mounting solution or just with a base for parking an existing mount. Fixation screws and dowels not included.

M3D Wall Mount 11.5000

Wall mount to operate a device that is attached to a wall. Includes two tubes (one of them L-shaped), a joint and an universal device socket (UDS). Includes slip out safeguards on all tubes and QuickShift levers on all joints. Requires mindful selection of a sturdy wall and appropriate dowels/screws (not included).



Wall mounts



M3D Wall mount Base Perpendicular

11.1051

Wall mounting plate with base for mounting Monty 3D tubes perpendicularly to the wall. Suitable for Monty 3D tubes. Fixation screws and dowels not included.



M3D Wall mount Base Parallel

11.1052


Wall mounting plate with base for mounting Monty 3D tubes parallel to the wall. The QuickShift lever allows for a quick adjustment and position change. Fixation screws and dowels not included.



M3D Storage Wall mount


11.5001

This wall mount unit consists of a wall mounting plate with base and a sturdy M3D tube. The mounting bracket allows Monty 3D wheelchair mounts to be removed from the wheelchair and packed away on the wall when their usage on the wheelchair is not needed. Fixation screws and dowels not included.



Parallel

to the wall



Perpendicular

to the wall



Socket

as mount storage